



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Surgical suction NOUVAG Vacuson 40 (New)

New

Swiss production

The advantages of the VACUSON 40 Sucker:

Economical due to the excellent price-performance ratio,

Precise and smoothly adjustable suction capacity,

Solid, reliable and proven technology,

Low noise and vibration levels,

Oil-free piston pump,

The Vacuson 40 device is designed for aspiration and removal of fluids and intraoperative contents with the possibility of adjusting the suction power, which allows adjusting the operation of the device to the needs of the procedure,

The mammal has a pressure gauge to control the settings,

The piston technology used in the Vacuson 40 mammals allows obtaining high vacuum parameters with a pressure of up to 0.9 bar and for long-term operation without any drop in parameters,

Thanks to their advantages, the devices can be used in both small and very large treatments,

Vacuson 40 has been equipped with an automatic anti-overflow system,

Standard equipment:

The VACUSON 40 control unit (0-1 bar manometer, vacuum regulator, 2 ring holders, 3 sockets: ON / OFF foot switch, foot controller and power cord, suction pipe with a set of connecting wires and a single bacterial filter),

Foot On / Off button,

Single-use bacterial filter (10pcs),

Power cord (5m),

Additional equipment (extra charge):

Jars of different capacity,

Trolley 'STAND' (2 shelves)

SANICAR mammal trolley (3 shelves),

Foot switch (vacuum regulator),

Technical data:

Air flow up to 40 l / min,

Vacuum: - 0.9 bar at 686 mmHG,

Energy consumption 185 W,

115/230 V, 50-60 Hz power supply,

Dimensions 290 x 280 x 360 mm,

Net weight 8.5 kg,

Has a current Technical Passport issued,

12 months warranty

