



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Surgical suction NOUVAG Vacuson 60 (New)**

New

Swiss production

The advantages of the VACUSON 60:

Economical due to the excellent price-performance ratio,

Precise and smoothly adjustable suction capacity,

Solid, reliable and proven technology,

Low noise and vibration levels,

Double piston pump without oil,

The Vacuson 60 suction pump motor drives two pistons to create a vacuum to achieve a capacity of 60 liters per minute,

The Vacuson 60 device is designed for aspiration and removal of fluids and intraoperative contents with the possibility of adjusting the suction power, which allows adjusting the operation of the device to the needs of the procedure,

The mammal has a pressure gauge to control the settings,

The piston technology used in the Vacuson 60 mammals allows obtaining high vacuum parameters with a pressure of up to 0.9 bar and for long-term operation without a drop in parameters,

Thanks to their advantages, the devices can be used in both small and very large treatments,

Vacuson 40 has been equipped with an automatic anti-overflow system,

Standard equipment:

VACUSON 60 control unit (0-1 bar manometer, vacuum regulator, 2 ring holders, 2 sockets: for ON / OFF foot switch, power cord, suction pipe with a set of connecting wires and a single bacterial filter),

Foot On / Off button,

Single-use bacterial filter (10pcs),

Power cord (5m),

Additional equipment (extra charge):

Jars of different capacity,

Trolley 'STAND' (2 shelves)

SANICAR mammal trolley (3 shelves),

Technical data:

Air flow up to 60 l / min,

Vacuum: - 0.9 bar at 686 mmHG,

Power consumption 250 W,

115/230 V, 50-60 Hz power supply,

Dimensions 290 x 280 x 360 mm,

Net weight 12 kg,

Has a current Technical Passport issued,

12 months warranty

