



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



RODENSTOCK R 7 vision tester (Reconditioned)

Used, very good technical condition,
German production,
Power supply: 230 V,
Frequency: 50-60 Hz,
Power: 15 W,
Intended, among others for tests of machine operators, electricians, etc.
Adapted to all kinds of spectacle frames,
Adjusting the angle of inclination,
Adjustable aperture left / right eye, binocular,
5 types of stereoscopic tests,
Eyepiece: 0, +1.0 D
Distance adjustment: 40cm, 55cm,
Dimensions of the device: 18 x 39 x 30/45 cm,
Has a current Technical Passport issued,
6-month warranty

