



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical monitor BEACON C61W + (New)

New,
Integrated, dual screen 6MP,
Monitor features:
DISPLAY - Built-in multi-image display modes,
DICOM - Compliance with the DICOM standard,
Application:
Radiografia Komputerowa / Direct Radiography - CR / DR,
Computed tomography - CT,
Magnetic resonance imaging - MRI,
Pozytonowa tomografia emisyjna / Tomography - PET / CT,
X-RAY X-ray,
Digital subtraction angiography - DSA,
Diffusion spectrum imaging - DSI,
3D medical images,
Nuclear medicine,
CAD medical software,
Reconstruction of images,
Technical data:
Size: 31.1",
Color depth: 1.07 billion,
Resolution: 3280 x 2048,
Brightness: 800 cd / m²,
Contrast: 1000: 1,
Viewing angle: 170° pion / 170° level,
Signal inputs:
DVI-D: 24 pin,
DisplayPort: 20 pin,
With table base,
Certificates: EN 60601, CCC, FCC, CE, UL,
Has a current Technical Passport issued,
24-month warranty

