



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BEACON S271T medical monitor (New)

New,

A color medical monitor for surgery and endoscopy,

Monitor features:

MULTI - Built-in multi-image display modes,

COLOR - Perfect color reproduction,

PERSPEKTIVE - Wide viewing angles,

BSS - System for stabilizing the matrix backlight,

Application:

Endoscopy,

Operating rooms,

Intensive care units,

Technical data:

Size: 27",

Color depth: 1.07 billion,

Resolution: 1920 x 1080

Brightness: 700 cd / m²,

Contrast: 1000: 1,

Viewing angles: 178° pion / 178° level,

Signal inputs:

DVI-D: 24 pin,

HD-SDI: BNC,

YPbPr: BNC * 3,

RGBS: BNC * 4,

CVBS: BNC,

S-VIDEO: DIN4,

VGA: D-Sub 15 pin,

Optical fiber: FC (optional),

Signal outputs:

DVI-D: 24 pin,

HD-SDI: BNC,

YPbPr: BNC * 3,

RGBS: BNC * 4,

CVBS: BNC,

S-VIDEO: DIN4,

Capacitive, 10-point touch pad,

Cap on the dies: Professional protective glass

Certifications EN 60601, CCC, FCC, CE, UL,

Has a current Technical Passport issued,

24-month warranty

