



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BEACON M982T medical monitor (New)

New,
8M color medical monitor for consultation,
Monitor features:
LARGE FORMAT - Large format monitor screen,
BSS - System for stabilizing the matrix backlight,
DICOM - Compliance with the DICOM standard,
TOUCH - Easy to use touch screen,
PROTECTION - Professional protective glass,
INTERACTIVITY - Interactive digital whiteboard,
INTEGRATED - Integrated mobile design,
Application:
Comprehensive consultations,
Doctor's trainings,
Demonstrations of operations,
Remote consultation,
Technical data:
Size: 98 ',
Color depth: 1.07 billion,
Resolution: 3840 x 2160,
Pixel size: 0.4845 x 0.4845 mm,
Brightness: 500 cd / m ²,
Contrast: 1300: 1,
Viewing angles: 178 ° pion / 178 ° level,
Response time: 12 ms,
Signal inputs:
DisplayPort: 20 pin * 4,
DVI-D: 24 pins,
Own weight: 160 kg,
Power supply: AC 230 V,
Power consumption: maximum 500 W,
Energy saving: EPA, VESA,
MTBF: > 50,000 hours,
Support for touch screen and digital whiteboard,
Built-in two independent working stands,
Has a current Technical Passport issued,
24-month warranty

