



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Rehabilitation Batriatric Therapy table 3 Section Couches (ST3673-SEERSMEDICAL) (New)

New,  
English production,  
The height of the table and the backrest is adjusted by means of an electric actuator,  
The table is equipped with large mobile wheels (central locking with steering function) (LMWD),  
Frame construction designed for easy use of a mobile lift,  
The table as a standard has a hand-held remote control,  
Cylinders used in the TiMotion table,  
Technical data:  
Total length of the table bed: 191 cm,  
The width of the bed is 80 cm,  
Safe working load and lifting capacity (from minimum height) 325 kg (715 lbs)  
Adjustment of the back section in the range from 0 ° to 75 °,  
Inclination angle of the footrest from -20 ° to + 75 °,  
Table height adjustment range from 45 cm to 98 cm,  
Handheld remote control and face hole with a standard plug,  
Exclusive 80 cm upholstery as standard,  
Additional accessories (extra charge):  
Triple foot switch (ST6209)  
Has a current Technical Passport issued,  
Warranty 24 - months for electric actuators



TAPICERKI STANDARD + (ZA DOPŁATĄ)



TAPICERKI STANDARD (W CENIE)

