



## Diagnostic and treatment table Medicare Mobile Treatment 2 (SM2450-SEERSMEDICAL) (New)

New,  
English production,  
Table height adjusted by means of a hydraulic cylinder,  
Adjustment of the backrest position by hand using a gas spring,  
The table is equipped with large individually braked wheels (LIBC),  
Frame construction designed for easy use of a mobile lift,  
The standard table has been equipped with a paper holder (6005) and side railing (6002),  
Frame coated with an antimicrobial powder coating,  
Cylinders used in the TiMotion table,  
Technical data:  
Total length of the table's bed: 194 cm,  
The width of the bed is 63 cm made of foam,  
Safe working load and lifting capacity 240kg (528 lbs),  
Backrest adjustment in the range from -25 ° to 75 °,  
Table height adjustment range from 46 cm to 99 cm,  
Large Tente® Ø125mm wheels,  
Additional accessories (extra charge):  
Upholstered side barriers - pair (6002P),  
Face hole with a plug in the headrest (6004),  
Upper limb support - Armrest (6006),  
Removable headboard - Cushion (6007),  
Base protection cover for long base (6053),  
Drip hanger (6054),  
Has a current Technical Passport issued,  
Warranty 24 - months for hydraulic cylinders and gas springs

