



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Trolley for transporting and washing the patient Shower Trolley (ST7800-SEERSMEDICAL) (New)**

New,  
English production,  
Table height regulated electrically from the mounted battery,  
The trolley has large wheels with a central lock that ensure high maneuverability,  
As standard, the stroller is padded with a PVC-coated foam with side props,  
The side barriers fold to make it easier to carry the patient,  
Easy drainage thanks to the outlet connected to the bottom of the base, connected to the drain pipe,  
On all corners of the trolley bumpers are mounted that provide protection during transport (bumpers can be removed in this case the width of the trolley is 65 cm),  
The height of the trolley can be lowered to the height of the wheelchair, which makes it easier to translate the patient,  
Frame construction designed for easy use of a mobile lift,  
Frame coated with an antimicrobial powder coating,  
Cylinders used in the TiMotion table,  
Technical data:  
Total length of the table bed: 195 cm,  
The width of the trolley with buffers 74 cm,  
Safe working load and lifting capacity 240 kg (528lbs),  
Height adjustment range from 41 cm to 94 cm with a manual switch,  
The SM 7800 with battery contains a spare battery and charging unit,  
Additional accessories (extra charge):  
Spare battery (6092)  
Replaceable PVC cladding (6093),  
Exchangeable PVC cushion (6094),  
Has a current Technical Passport issued,  
Warranty 24 - months for hydraulic cylinders

