



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Injection needle MEDWORK INJ1-A1-07-5-23-230 (Ex-Display)

New, exhibition,
German production,
Diameter: 0.7mm,
Length: 230 cm,
Operating channel: 2.3 mm,
One-time,
Expiration date: 04.2018

