



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DA803 dental compressor (New)

New,
Power supply: 230 V,
Frequency: 50 - 60 Hz,
Power: 1.5 kW,
Oil-free,
Silent diaphragm compressor with advanced microcomputer control:
Monitoring of work status in real time,
Monitoring of voltage, current, pressure, etc. in real time,
Switching and stopping pressure regulated via a digital display,
Intelligent fault warning system,
Independent control of each engine to ensure the continuity of work on failure, one of the engines,
Air flow: 304 l / min,
Capacity at 7 Bar: 100 l / min.,
Tank capacity: 50 liters,
Sound intensity: less than 70 dB,
Weight: 54 kg,
Overall dimensions: 675 x 480 x 675 mm,
Has a current Technical Passport issued,
12-month warranty

