



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



HEWLETT PACKARD M3000A module (Reconditioned)

Used, very good technical condition,
German production,
The module has:
ECG / Resp,
SpO2
NBP
PRESS / Temp
Compatible with:
AGILENT Model M3,
PHILIPS Model INTELLI VUE,
HEWELETT PACKARD VIRIDIA M3,
Overall dimensions: 10 x 19 x 5.5 cm,
Has a current Technical Passport issued.
1-month warranty

