



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Cryosurgery device - portable Freezpen CL-FP-18 (New)**

New,  
Belgian production,  
Freezpen is manufactured by Clinic6,  
It is an innovative device that increases the overall improvement in the quality of treatment through the direct use of Nitrogen Oxide at  $-89^{\circ}\text{C}$ , penetrates in 1 mm for 5 seconds to a maximum depth of 5 mm,  
Ease of use is an incentive for use by family doctors, aesthetic medicine and pediatricians,  
Freezpen offers stress-free treatment with quickly visible improvement effects,  
High precision Freezpen limits damage to healthy skin, minimizing discomfort in patients,  
Effective result obtained thanks to the first treatments, thanks to the high pressure of Freezpen, achieving deeper penetration,  
The device has a calming feather shape, the advantage is optimal safety with minimal discomfort and no side effects,  
Very economical for cheap cryogenic capsules that can heal up to 12 patients,  
The camera uses cartridges weighing 16 grams,  
The device is the optimal solution for treatment:  
warts,  
condyle,  
hemangiomas,  
skin markers,  
cornea and all other mild skin disorders.  
The camera includes:  
3 applicators:  
1 applicator with blue coding (for the treatment of lesions of diameter 1-3 mm)  
1 applicator with white coding (for treatment of lesions diameter 2-6 mm),  
1 applicator with green coding (for the treatment of lesions diameter 4-8 mm),  
Instructions for use,  
Tweezers for filter placement.  
Has a current Technical Passport issued,  
12 months warranty

