



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **WEARCARE WT1 wireless bluetooth thermometer (New)**

New,  
Production P.R.C,  
Power supply: 3V 225 mAh battery,  
Constant monitoring of body temperature in real time,  
Alarm when the set temperature zone is exceeded and Bluetooth connection is interrupted,  
In addition to high temperature alarms, the Temp Sitter function can also display temperature and trend graphs in real time using a graphic grid and record treatments such as antipyretic or physical cooling,  
Bluetooth Intelligent data transfer to a smartphone,  
You do not need to charge or turn on / off the device, very low power consumption, you can use it for 1 year without changing the battery (8 hours / day),  
Application software: BT4.0 and IOS 6.0 above, Android versions 4.3 and higher,  
You can choose 1 of 4 colors,  
Included:  
Thermometer + case,  
10 self-adhesive pads,  
Leash,  
User manual in English,  
Weight: 9.6 g,  
Overall dimensions: 51 x 31.8 x 8.5 mm,  
6-month warranty

