



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Medical couch ENGO S01 (New)**

New,  
Of Polish origin,  
ENGO S01 couch, foldable to the wall, is made of powder-coated steel tubes,  
The couch is equipped with a fixed height and one-part sleep surface with a face hole and a plug,  
The ENGO S01 couch is intended for rehabilitation centres, hospitals, health centres and rehabilitation, treatment and doctor's surgeries,  
It is possible to order upholstery in a colour of your choice free of charge, except for "PRESTIGE" colours which are available at an extra charge,  
As additional options are available (at extra charge):  
Hanger for the sleeper,  
Technical parameters:  
Length 190 cm,  
Width 55 or 70 cm,  
Height 55 or 70 cm,  
Weight 35 kg,  
Maximum load 150 kg,  
It has a current Technical Passport,  
Warranty: 24 months,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE

