



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Medical couch ENGO K (New)**

New,

Made in Poland,

It has a simple and stable frame made of powder-coated steel tubes,

The ENGO K couch is of fixed height and has a two-section sleep surface with adjustable headrest with a rastomat mechanism,

The standard version has a cut-out on the left side to facilitate cardiological examination,

It is possible to order upholstery in a colour of your choice free of charge, except for "PRESTIGE" colours which are available at an extra charge,

Additional options are available (at extra charge):

Hanger for the sleeper,

Castors with brake,

Technical parameters:

Length 190 cm,

Width 55 or 70 cm,

Height 55 or 70 cm,

Weight 35 kg,

Maximum load 150 kg,

Headrest angle adjustment 0° to +45°,

Side angle adjustment: 0° to 90°,

It has a valid Technical Passport,

Construction warranty 24 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

