



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Defibrillator PROGETTI Rescue Life (7 'Version) (New)

New,
Italian production,
Portable (with handle),
Progetti Rescue Life is a professional manual clinical defibrillator intended for advanced resuscitation procedures by trained personnel,
Software: multilingual voice and text prompts,
Memory: Internal flash memory 1 Gb, up to 300 hours of ECG, events, review recording,
Menu and operation in English,
Technical data
Screen:
Large 7" color TFT display,
ECG signal: from spoons or 3-lead cable: I, II, III (option 6- or 12-lead leads),
SpO2 and trace,
Heart rate,
ECG Monitor:
Sweep: 5, 12.5, 25, 50 mm/sec,
Amplitude (feature) 5, 10, 20 mm/mV,
Heart rate - 20 to 300 bpm ($\pm 2\%$),
Alarm: ability to set HR alarm in range 20-250 bpm,
Printer:
integrated thermal printer for ECG traces and events, printout with HR/SpO2 values; paper 58mm.
Power:
High-capacity NiMH battery with short charging time (95% within 2h),
Mains/charging: 100-230 V, $\pm 10\%$ 50/60 Hz,
Standard equipment:
Rescue Life 7 defibrillator with NiMH battery and spoons,
Bag,
3-core cable,
1 roll of paper,
Optional equipment/parameters (for an extra charge):
SPO2 Module,
NIBP Module,
External pacemaker module (pacing rate 30-180/min; amplitude 0-180mA, time 22ms),
Pediatric spoon set,
Car charger,
Disposable electrodes,
10-core cable (12-lead),
Pulse type:
BTE (two-phase exponential pulse); impedance compensated,
Energy: Manual: 1-230J,
Energy: AED: 150J - constant,
HiCAP technology (Large Storage Capacitor),
Defibrillation energy selection: 1 - 230J - multifunction knob on panel,
Charging time:
100 discharges with charged batteries,
Synchronisation mode of discharge with R (cardioversion),
Charging/defibrillation with a button on the spoons,
Memory:
Internal flash memory 1 Gb, up to 300 hours of ECG, events, review recording,
Dimensions:
369mm x 240mm x 340 mm,
Weight:
5.5 kg (with battery),
Includes current Technical Passport,
2-year warranty,
Financing options (Poland only): Installments, Leasing, Loan,

