



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **AED Training Defibrillator PROGETTI Rescue SAM PRO (New)**

New,  
Italian production,  
Portable (equipped with a handle),  
The two-phase,  
RESCUE SAM PRO is a semi-automatic, external training defibrillator designed to provide a simple, uncomplicated service, voice prompts inform about the next steps during the resuscitation process,  
Technical data:  
OPERATING MODE: semi-automatic,  
PROTOCOL: Voice prompts and indicators guide you through the defibrillation process,  
VOICE INSTRUCTIONS: detailed voice messages from the guide during the whole defibrillation process,  
CONTROLS:  
shock button,  
information button,  
on / off button  
INDICATORS: LED (various colors), information button,  
SCENARIOS - 5 preset:  
1 Flickering of chambers with single shock conversion,  
2 Flickering of chambers with multiple shocks required for conversion,  
3 Flickering of chambers with two shocks necessary for conversion,  
4 Non-defibrillation rhythm,  
5 Vibrations of chambers with three shocks necessary for conversion,  
TRAINING ELEKTRODES: 3 pairs of electrodes (1 pair for adults and 2 pairs for children - electrodes are reusable),  
BATTERY: 9V NiMH battery for long-lasting working time,  
STATIC ELECTRICITY: <5uA,  
REMOTE CONTROL: DC 3V / 2 AAA battery,  
DIMENSIONS: 290x280x90 mm,  
WEIGHT: 2 kg,  
Included:  
defibrillator,  
Remote,  
Training electrodes and adapters,  
Charger,  
User manual,  
Bag for electrodes,  
Additional equipment (extra charge):  
Bag,  
Has a current Technical Passport issued,  
24-month warranty

