



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



AED Training Defibrillator PROGETTI Rescue SAM PRO (New)

New,
Italian production,
Portable (equipped with a handle),
The two-phase,
RESCUE SAM PRO is a semi-automatic, external training defibrillator designed to provide a simple, uncomplicated service, voice prompts inform about the next steps during the resuscitation process,
Technical data:
OPERATING MODE: semi-automatic,
PROTOCOL: Voice prompts and indicators guide you through the defibrillation process,
VOICE INSTRUCTIONS: detailed voice messages from the guide during the whole defibrillation process,
CONTROLS:
shock button,
information button,
on / off button
INDICATORS: LED (various colors), information button,
SCENARIOS - 5 preset:
1 Flickering of chambers with single shock conversion,
2 Flickering of chambers with multiple shocks required for conversion,
3 Flickering of chambers with two shocks necessary for conversion,
4 Non-defibrillation rhythm,
5 Vibrations of chambers with three shocks necessary for conversion,
TRAINING ELEKTRODES: 3 pairs of electrodes (1 pair for adults and 2 pairs for children - electrodes are reusable),
BATTERY: 9V NiMH battery for long-lasting working time,
STATIC ELECTRICITY: <5uA,
REMOTE CONTROL: DC 3V / 2 AAA battery,
DIMENSIONS: 290x280x90 mm,
WEIGHT: 2 kg,
Included:
defibrillator,
Remote,
Training electrodes and adapters,
Charger,
User manual,
Bag for electrodes,
Additional equipment (extra charge):
Bag,
Has a current Technical Passport issued,
24-month warranty

