



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



UNIKA PLUS dryer (New)

New
Italian production,
Tecno-Gaz has created a unique and extremely rational mechanical device for drying sterilized materials,
This device allows for obligatory drying of instruments and avoids the use of instruments (especially sharp instruments),
The dryer guarantees perfect drying, thus preventing possible corrosion and white marks on the instruments,
UNIKA is an indispensable, very useful device for preparing instruments. The average time required for good drying is only about four minutes.
Advantages of the dryer:
Completely made of stainless steel,
Equipped with a double convector for transferring hot air to the drying chamber,
Special tray supplied with the device,
Timer,
Technical data:
Overall dimensions:
Width: 200 mm,
Height: 200 mm,
Depth: 305 mm,
Weight: 4.5 kg,
Power supply: 220 V - 50/60 Hz,
Power consumption: 1 kW,
Has a current Technical Passport issued,
12 months warranty

