



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Manual welder ONE (New)

New,  
Italian production,  
ONE is a new generation sealing machine that has been developed to reduce envelope packing time, monitor welding quality and seal all types of paper for sterilisation, including folded paper,

Advantages of the sealer:

12 mm sealing strip,  
32 cm sealing strip,  
Simple paper insertion,  
Active electronic control of several heating elements,  
Visible cutting area,  
Warning of incorrect sealing,  
Can be placed anywhere,

Technical data:

Dimensions:

Width: 460 mm,  
Depth: 400 mm,  
Height: 170 mm,  
Weight: 8 kg,  
Power consumption: 320 W,  
Power supply: 230 V~ - 50 / 60 Hz,  
Seal width: 12 mm,  
Sealing length: 320 mm,  
With issued valid technical passport,  
Warranty: 12 months,

Finance options (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

