



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Dea Lux whitening lamp (New)**

New

Italian production,

Dea Lux uses a wavelength from 480 to 520 nm,

Additionally, it has a special screen that allows absorption of ultraviolet and infrared light, thus preventing the heating of the pulp and mucous membranes and avoiding / reducing their sensitivity,

A special optical fiber screen controls the light and ensures the concentration of rays, which allows maximum efficiency and shortens the treatment time,

The large liquid crystal display enables the user to monitor all functions of the device and display each one,

Dea Lux takes up very little space,

It has a solid, stable, balanced base with wheels, which makes it easy to move around. The full rotation system allows for any positioning.

Technical data:

Dimensions (X, Y, Z): 1000 x 390 x 315 mm,

Weight: 15 kg,

Power supply: 230 V, 50 Hz,

Number of phases: 1,

Power consumption: 50 W,

Radiation wavelength: 495 nm,

Output power: 30 W,

Has a current Technical Passport issued,

12 months warranty

