



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



VITAL TEST pulse oximeter (New)

New,
Italian production,
Vital Test is a professional pulse oximeter for use in dentistry,
Simple, with many functions, it can be equipped with a special holder for mounting on the device, wall or stand. This option allows continuous control of the instrument, making it convenient for continuous use,
The advantages of the pulse oximeter:
Simple and easy to use,
Large LCD display with backlight function,
The ability to set and turn off alarms for lower and upper thresholds of saturation and heartbeat,
Displaying trends,
Visualization of the pletic wave,
Memory of obtained measurements,
Sensor fault indication or disconnection,
Additional equipment (extra charge):
Holder for mounting the pulse oximeter at the wall,
The Vital Test can be delivered with a programmable interface for connection to a computer, for recording and printing of all measured parameters,
Technical data:
Dimensions (LxWxH): 160 x 70 x 376 mm,
Weight: 165 (g) (without battery),
Battery-powered: 4 x 1.5 V ~ LR6 AA (48 hours),
Has issued the current Technical Passport:
12 months warranty

