



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Autoclave series YESON E-8I (New)

New,

A modern Yeson E-series autoclave with a capacity of 8 liters,

The initial fractional vacuum (class B - medical) allows you to sterilize a variety of unwrapped and packaged tools - solid, porous, narrow-clearance articles in accordance with the PN-EN 13060-2015 classification,

Intended for:

doctors' offices,

dental,

cosmetic,

veterinary,

clinics,

hospitals,

Features and benefits of the E-series autoclave:

Chamber made of stainless steel,

Class B - fractional vacuum,

Built-in printer,

USB connector for pendrive,

LCD display,

Touch control buttons,

Polish and English menu,

Block against accidental opening,

Automatic fuses,

9 work programs,

2 test programs,

Delayed start function,

Pre-heating function,

Password protection function

Standard equipment:

Enlarged basket with 3 trays,

Power cord,

Holder for trays,

Shutter cable,

Polish instruction,

Technical data:

Capacity: 8l,

Power 1750 W,

Maximum temperature 137 ° C,

Dimensions of the housing (height x width x depth): 320 x 365 x 535 mm,

Chamber dimensions Ø 170 x 320 mm,

Weight: 50 kg,

European CE certificate,

Has a current Technical Passport issued,

Warranty 24 months

