



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Autoclave YESON E series 22l (New)

New,

Modern Yeson E-series autoclave with a capacity of 22 litres,

Pre-vacuum fractionated (class B - medical) allows to sterilize a variety of unpackaged and packaged instruments - solid products, porous products, objects with narrow clearance according to classification PN-EN 13060-2015,

Designed for:

doctor's surgeries,

doctor's offices, dentist's

cosmetics,

veterinary surgeries,

clinics,

hospitals,

Features and benefits of E-series autoclave:

Stainless steel chamber,

Class B fractionated vacuum,

Built-in printer,

USB connector for pen drive,

LCD display,

Touch control buttons,

Polish and English menu,

Anti-accidental opening lock,

Automatic fuses,

9 operating programmes,

2 test programmes,

Delayed start function,

Pre-heating function,

Password protection function

Standard equipment:

Enlarged basket with 3 trays,

Power cord,

Tray holder,

Drain hose,

Polish manual,

English instructions,

Technical data:

Capacities: 22l,

Power: 2.3 kW or 2.4 Kw respectively,

Maximum temperature 137° C,

Unit dimensions: 510mm x 680mm x 430mm,

Dimensions of the chamber (22l) Φ 250 x 450 mm,

Weight: 54 kg,

European CE Certificate,

It has a valid technical passport,

Warranty: 24 months

