



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Solvent recycler water Steel 4 I (New)**

New

4L water distiller with a stainless steel interior,

All elements of the water cycle (tank, head interior, coil, filter casing, jug cap) made of SU304 steel to meet the expectations of the most demanding customers,

Operation based on a cooler with forced air circulation,

Application:

For the production of distilled water of higher quality than chemically demineralised, in order to ensure reliable operation of various devices,

Irreplaceable in dental offices, laboratories,

The obtained water is distilled to:

Autoclaves and dental units,

Air humidifiers,

Dilution of alcohols and other liquids,

Standard equipment:

Destylarka,

Water tank with connector and filter,

Power cord,

Cleaning powder,

User manual in Polish,

Technical data:

230V (220-240V) power supply,

Power consumption 750W,

Water tank capacity 4 liters,

Capacity - about 1 liter per hour,

Has a current Technical Passport issued,

Warranty 12 months

