



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Seal 320 impulse sealer (New)**

### **New**

A modern impulse welding machine with reduced energy consumption, equipped with an electronically controlled heating element. Maximum length of the weld 320mm,

It allows the preparation of paper and foil packages intended for sterilization of tools in an autoclave,

Conducting sterilization in packages is considered a particularly effective, safe and convenient method. It allows you to store sterile tools for a long time after the process,

Features and advantages of the welding machine:

Welding machine can be used to close ready-made packaging or sleeves in rolls with a maximum total width of 320 mm,

Precise cutting off the necessary section of the sleeve allows a two-sided knife installed in a safe way for the operator. Reinforced blades do not need to be replaced even during intense operation,

The housing with a modern look from solid, easy to clean plastic and stainless steel,

LED device status indication,

Automatic transition into standby mode,

Welding with a width of 12 mm guarantees tightness of packets,

Even pressure over the entire length of the weld,

Technical data:

Voltage 230 V / 50 Hz,

100 W power,

Maximum length of the weld 320 mm,

The width of the weld 12 mm,

Heating time 3 minutes,

Dimensions: D. 460 x Sz. 400 x W. 170 mm,

The weight of the device is 7 kg,

Has a current Technical Passport issued,

12 months warranty

