



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Welding machine F-300 YESON (New)**

New

Features and advantages of the welding machine:

The device is helpful in the preparation of paper and foil packages intended for sterilizing tools in an autoclave,

Conducting sterilization in packages is considered a particularly effective, safe and convenient method,

The welder can be used to close ready-made packages or sleeves in rolls of different widths,

Precise cutting off the necessary section of the sleeve is possible with a double-blade knife,

Temperature control of the heating element helps to obtain very good quality welds,

The robust stainless steel housing and the metal construction of the device guarantee stability and durability,

The pressure of the sleeve to the heating element made of heat-resistant silicone,

Technical data:

Voltage 230 V / 50 Hz,

500 W power,

Maximum length of the weld 250 mm,

The width of the weld 12 mm,

The weight of the device is 5 kg,

Has a current Technical Passport issued,

12 months warranty

