



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Segment - leg rest ESCHMANN TA-0301001 (Reconditioned)

Used, perfect technical condition,
English production,
Very light construction,
Mounting distance: 48.5 cm,
Maximum load: 40 kg,
Overall dimensions: 56 x 86 x 10 cm,
Has a current Technical Passport issued,
1-month warranty

