



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Universal headrest MAQUET (Reconditioned)

Used, very good technical condition,
German production,
Manual adjustment in various positions,
Manual adjustment of the headrest spacing,
Dimensions of the headrest: 23/30 x 20.5 x 5 cm,
Has a current Technical Passport issued,
1-month warranty

