



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Hawkmед HK-300 pump (New)**

New

Features and benefits of the pump:

Compact and lightweight design saves space and is convenient when moving the patient,

The panel lock function helps to prevent unauthorized changes to the settings of each instrument,

Large and colorful LCD display,

History of 2000 infusions,

RS232 interface (optional),

Real-time display,

Adjustable buzzer volume (3 levels),

Technical data:

Flow rate range: 1 ~ 400 ml / h,

The volume to be administered (VTBI): 0 ~ 9999 ml,

Infusion range ( $\Sigma$ ): 0 ~ 36000 ml,

Accuracy of the infusion:  $\pm 10\%$ ,

The right food bag: Serve various brands of nutritional preparations,

Bolus flow rate: 400 ml / h,

Occlusion pressure: 3 adjustable occlusion pressure settings: low, medium and high,

Alarms: Visual and audible alarms: Open door, occlusion, end of infusion, empty, start-up reminder function, discharged battery, discharged battery, failure, etc.

Power supply: AC: 100 ~ 240V, 50 / 60Hz DC: 12V  $\pm$  1V,

Battery: Lithium-polymer battery, 7.4V, 1900mAh. It can work for more than 6 hours when fully charged at 25 ml / h,

Operating mode: continuous,

Dimensions: 145 x 100 x 120 mm (length x width x height),

Weight:  $\leq 1.4$ kg,

Power consumption: 25VA,

Classification: Class I, type CF,

Water tightness - IPX3 standard,

Working conditions: Temperature: 5 ° ~ 40 °, relative humidity: 10 ~ 95%, atmospheric pressure: 86 ~ 106 kPa,

Transportation conditions: Temperature -20 ° ~ + 60 °, relative humidity: 10 ~ 95%, atmospheric pressure: 50 ~ 106 kPa,

Storage conditions: Temperature -20 ° ~ + 45 °, relative humidity: 10 ~ 95%, atmospheric pressure: 50 ~ 106 kPa,

Has a current Technical Passport issued,

Warranty 12 months

