



Laboratory centrifuge BIO RAD CONTRIFUGE 24 (Reconditioned)

Used, very good technical condition,

German production,

Power supply: 230 V,

Frequency: 50-60 Hz,

Power: 70 W,

BIO RAD Centrifuge 24 is a stationary laboratory centrifuge with a microprocessor control designed for testing red blood cells agglutination in gel identification cards,

Centrifuge designed for 24 ID cards,

Centrifugation of preparations at a constant speed: 900 rpm,

Constant working time of the centrifuge: 10 min,

It has a brake,

It has an electromagnetic lock for opening the cover,

Built-in automatic balance adjustment,

It is characterized by a very quiet mode of operation,

Included:

The angle rotor has 24 slots for plates (cards) with a width of 5 mm and a length of 70 mm,

The diameter of the chamber: 32 cm,

Overall dimensions: 37 x 40 x 23 cm,

Has a current Technical Passport issued,

1-month warranty,

Kat 02

