



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Light source OLYMPUS CLD-S (Reconditioned)

Used, very good technical condition,
Japanese production,

Power supply: 230 V,

Frequency: 50-60 Hz,

Adapted to rigid endoscopes,

Source of light:

150 W halogen meta light bulb - working lighting (bulb consumption at the level of 40%),

Color temperature: 5600 K,

5-step adjustment of light intensity,

Overall dimensions: 37 x 38 x 13 cm,

Has a current Technical Passport issued,

Warranty for 1 month

