



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Rehabilitation bed FAMED ZYWIEC LR-12 (Reconditioned)**

Used, perfect condition,

Polish production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power consumption: 300 W,

Rehabilitation bed LR-12 is intended for adult hospital patients during diagnosis, monitoring and treatment of diseases,

The autoregression function in the backrest segment reduces pressure on the patient's body in the lumbar and thigh regions,

Trendelenburg and anti-Trendelenburg tilt adjustment by gas spring, adjustment levers easily accessible from both sides of the bed,

Technical data:

Bed length: 2160 mm,

Bed width: 965 mm,

Length/width of mattress: 1960 x 865 mm,

Minimum lying height: 430 mm,

Maximum bed height: 805 mm,

Backrest segment lift angle: 0 - 70°,

Thigh segment elevation angle: 0 - 37°,

Shank segment elevation angle: 20°,

Trendelenburg position: 14°,

Anti-Trendelenburg position: 20°,

Backrest self-preservation function: 90 mm,

Permissible load: 170 kg,

The diameter of the wheels: 150 mm,

The CPR lever of the backrest segment, the quick release lever of the segment is used especially in emergency situations (e.g. resuscitation),

enables quick levelling,

Beds equipped with round bumpers,

Included:

Exhaust frame that includes:

The leg of the frame from the side of the patient's head: 2 pieces.,

Frame leg from the side of the patient's legs: 2 pcs.,

Leg jointing bar: 2 pcs,

Sliding bar: 2 pcs,

I-barrier: 1 pcs,

Crossbar II: 2 pcs,

Roller handle: 1 pc.,

Jaw chuck: 10 pcs.,

No mattress included - can be retrofitted for an additional fee,

Mobile on wheels with centralised brake and directional travel function,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan.

