



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Set of ophthalmic tools (version C) (Reconditioned)**

Used, very good technical condition,  
They are used to carry out microsurgical procedures and surgeries on the eye,

Set contains:

4 pieces of slotted knives, Double angled blade at an angle of 45 °,

2 pieces 1.2 mm,

1 piece 2.8 mm,

1 piece 2.65 mm,

1 piece of precision tools,

4 pieces of tweezers,

1 eyelid retractor,

Scissors,

Plastic casing for sterilization

Autoclavable plastic with silicone mat; the base, lid and mat are perforated to penetrate the steam,

External dimensions: 265 x 160 x 25 mm,

Has a current Technical Passport issued,

1-month warranty

