



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Operating light LED EMALED 200 D (Ceiling) (New)

New,
German production,
Ceiling mount lamp,
Advantages and characteristics of the lamp:
Works well in examination and treatment rooms,
Light with high contrast and low energy consumption,
shadowless lighting and good color reproduction,
Low power consumption of test light with a capacity of just 10-15 W,
Long lifetime of LEDs of minimum 50,000 hours ensures low operating costs,
Low service costs and a long service life ensure safe operation and low maintenance costs,
Light can be set in any position thanks to spring arm with adjustable height and 360 ° rotation,
high mobility and good lamp body position,
natural color reproduction,
technical data:
maximum light intensity (per 1 m): 50,000 lux,
lighting intensity: 30 - 100%,
Operating distance: 70 - 140 cm,
Color temperature: 4 500 K,
Coefficient of ret color rendering 95 Ra,
Durability of LEDs 50,000 h,
Diameter 24 cm,
7 LEDs,
Fixed focal length,
Main voltage: 230 V,
Frequency 50 - 60 Hz,
10VA lamp head power,
Header weight lamps 1 kg,
Total weight 12 kg,
Included power supply and a set for ceiling mounting,
680mm suspension tube (can be ordered with a 1500mm suspension pipe for an additional fee),
Has a current Technical Passport issued,
Warranty 24 months

