



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **Operating light LED EMALED 200 D (Ceiling) (New)**

New,  
German production,  
Ceiling mount lamp,  
Advantages and characteristics of the lamp:  
Works well in examination and treatment rooms,  
Light with high contrast and low energy consumption,  
shadowless lighting and good color reproduction,  
Low power consumption of test light with a capacity of just 10-15 W,  
Long lifetime of LEDs of minimum 50,000 hours ensures low operating costs,  
Low service costs and a long service life ensure safe operation and low maintenance costs,  
Light can be set in any position thanks to spring arm with adjustable height and 360 ° rotation,  
high mobility and good lamp body position,  
natural color reproduction,  
technical data:  
maximum light intensity (per 1 m): 50,000 lux,  
lighting intensity: 30 - 100%,  
Operating distance: 70 - 140 cm,  
Color temperature: 4 500 K,  
Coefficient of ret color rendering 95 Ra,  
Durability of LEDs 50,000 h,  
Diameter 24 cm,  
7 LEDs,  
Fixed focal length,  
Main voltage: 230 V,  
Frequency 50 - 60 Hz,  
10VA lamp head power,  
Header weight lamps 1 kg,  
Total weight 12 kg,  
Included power supply and a set for ceiling mounting,  
680mm suspension tube (can be ordered with a 1500mm suspension pipe for an additional fee),  
Has a current Technical Passport issued,  
Warranty 24 months

