



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Operating light LED EMALED 200 W (Wall) (New)**

New,  
German production,  
Lamp for wall mounting,  
Advantages and characteristics of the lamp:  
Works well in examination and treatment rooms,  
Light with high contrast and low energy consumption,  
No moving parts in the lamp head,  
Shadowless lighting and good color reproduction,  
Low power consumption of test light with a capacity of just 10-15 W,  
Long lifetime of LEDs of at least 50,000 hours ensures low operating costs,  
Low service costs and a long service life ensure safe operation and low maintenance costs,  
Aluminum housing covered with powder,  
High mobility and a good setting of the lamp body,  
Barely noticeable heating of the lamp head,  
Natural color reproduction,  
3-step brightness adjustment,  
Two auto clivable holders are included (additional holders can be ordered at any time),  
Specifications:  
Maximum light intensity (per 1 m): 50,000 lx,  
Lighting intensity: 30 - 100%,  
Operating distance: 70 - 140 cm,  
Color temperature: 4 500 K,  
Color reproduction ratio 95 Ra,  
Durability of LEDs 50,000 h,  
Diameter 24 cm,  
7 LEDs ,  
Fixed focal length,  
Main voltage: 230 V,  
Frequency 50 - 60 Hz,  
10VA lamp head power,  
Lamp head weight 1 kg,  
Total weight 10 kg,  
Included wall power supply with electrics, spring arm, boom ,and light head EMALED 200,  
Has a current Technical Passport issued,  
Warranty 24 months

