



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Operating light LED EMALED 200 F (Mobile) (New)**

New,  
German production,  
Mobile lamp (base with power supply),  
Advantages and characteristics of the lamp:  
Works well in examination and treatment rooms,  
High mobility and good body positioning,  
Light with high contrast and low wear energy,  
No moving parts in the lamp head,  
shadowless illumination and good color reproduction,  
Low power consumption of test light with a power of just 10-15 W,  
Long lifetime of LEDs of at least 50,000 hours ensure low operating costs,  
low costs service and long service life ensure safe operation and low maintenance costs,  
Aluminum powder coated housing,  
High mobility and good lamp body positioning,  
Barely noticeable heating of the lamp head,  
Natural color reproduction,  
3-step brightness adjustment,  
Two are included auto clavable holders (additional handles can be ordered at any time),  
Data tech:: Maximum light intensity (per 1 m): 50,000 lx,  
Light intensity: 30 - 100%,  
Operating distance: 70 - 140 cm,  
Color temperature: 4 500 K,  
Color reproduction ratio 95 Ra,  
LED lifetime 50,000 h,  
Diameter 24 cm,  
7 LEDs,  
Fixed focal length,  
Main voltage: 230 V,  
Frequency 50 - 60 Hz,  
10VA lamp head power,  
1 pc lamp head weight,  
Weight total 20 kg,  
Included mobile base with 230V power supply, rotating base, spring arm, and head EMALED 200  
Has a current Technical Passport issued,  
Warranty 24 months

