



Laboratory centrifuge KENDRO / SORVALL LEGEND T (Reconditioned)

Used, very good technical condition,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 1250 W
Built-in numeric displays,
Foiled control panel,
Speed control: 300-15000 rpm, with the rotor used max speed 5650 rpm,
Time switch in the range: 1 sec - 9 h 59 minutes, or permanent work,
Three-step adjustment of the acceleration force,
Electric lock cover lock,
Induction brushless motor guarantees long-term operation without maintenance,
Centrifuge chamber made of stainless steel,
Chamber diameter: 46 cm,
Swinging rotor 75006444 which has 2 slots for microplate blocks,
Sound level: 66 dB (A),
Dimensions: (WxDxH) 551 x 666 x 362 mm,
User manual in English (PDF),
Has a current Technical Passport issued,
1-month warranty

