



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Diagnostic and treatment table SK-06 (New)**

New,  
Made in Poland,  
It is characterised by stable steel construction and column lifting mechanism, which eliminates movement of the table top during height adjustment,  
It has a 4-part sleep surface,  
The adjustment of the height and particular parts of the top (except for the bolster) is carried out with the use of an electric actuator controlled with a manual remote control.  
FANOS 02 has a wide range of adjustment of the sleep surface to the lying position, armchair, Pivot, Trendelenburg and Anti-Trendelenburg positions or any combination of four sleep surface segments,  
It is possible to order upholstery in a colour of your choice free of charge, except for "PRESTIGE" colours which are available at an extra charge,  
Table adjustment is carried out with the use of  
manual remote control as standard - serves for 4 adjustments  
Foot remote - serves for 4 adjustments (for an extra charge),  
There can be 2 remote controls in the table but it cannot be controlled from the frame. Only two controls can be installed in the table at the same time  
As additional options are available (for an extra charge):  
Armrests,  
Paper hanger,  
Foot remote control (instead of manual remote control),  
Technical data:  
Length: 200 cm  
Selectable width: 65/70/80 cm  
Height: 59-99 cm  
Headrest angle adjustment: -70° to 45° (by gas spring),  
Adjustment of the angle of the chest section: -17° to 70° (with electric actuator),  
Lumbar section angle adjustment: -25° to 20° (with electric actuator),  
Footrest angle adjustment: -40° to 25° (by electric actuator),  
Trendelenburg position: 0° to 25° (by electric actuator),  
Anti-Trendelenburg position: -17° to 0° (with electric actuator),  
Weight: 113 kg  
Maximum load: 200 kg  
It has a valid Technical Passport,  
Warranty: 24 months for the structure,  
Warranty on hydraulic and electrical accessories 12 months,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE

