



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Gynecological chair SELENE (New)

New,
Polish made,
It has three independent electric actuators, controlled by a hand remote control,
The chair is equipped with a set of euro rails mounted on the sides of the seat, a set of knee pads and a stainless steel bowl placed under the seat,

It is easy to keep clean and hygienic due to the built-in construction elements,
It is possible to order upholstery in a colour of your choice free of charge, except for "PRESTIGE" colours which are available at an extra charge,

For an extra charge, the armchair is available with the options shown in picture 7:

Armrests (1),
Removable footrest (2),
Paper hanger (3),
IV hanger (4),
Injection support (5) - upholstery can be changed to chair colour at extra cost,
Foot attachment (6),

Technical data:

Length: 135 cm,

Width: 57 cm,

Height: 55 - 95 cm,

Backrest angle adjustment: from 0° to 74° (by electric actuator),

Trendelenburg position: 0° to 19° (by electric actuator),

Anti-Trendelenburg position: -6° to 0° (by electric actuator),

Weight: 90 kg,

Maximum load: 200 kg,

It has a valid Technical Passport,

24 months construction warranty,

Warranty on hydraulic accessories 12 months,

Possible financing (Poland only): Installments, Leasing, Loan

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

