



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



---

### **MAQUET 3101.0234 holding hand cuff (Reconditioned)**

---

Used, perfect technical condition,

German production,

The strap can be removed for washing,

Included:

Mounting bracket for Maquet 1002.24C0 rail,

Has a current Technical Passport issued,

1-month warranty

