



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Patient arm strap MAQUET 1002.24C0/01 (Reconditioned)

- Used, perfect technical condition,
- German production,
- The strap can be removed for washing,
- **Included in the set:**
 - Maquet 1002.24C0 rail mounting bracket,
- Has a valid Technical Passport issued,
- Guarantee:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Financing options (Poland only): Installments, Leasing, Loan,

