



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Nasal speculum 24mm pliers (Reconditioned)

Used, very good condition,
Polish production,
Made of stainless steel,
Intended for repeated use,
The possibility of sterilization,
Sight glass length 24 mm,
Model No. 0020

