



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BEACON C44W + medical monitor (New)

New,
Integrated, dual screen 4MP,
Monitor features:
Integrated two screens,
Built-in multi view mode,
System for stabilizing the matrix backlight,
Perfect color reproduction,
DICOM - Compliance with the DICOM standard,
Application:
Radiografia Komputerowa / Direct Radiography - CR / DR,
Computed tomography - CT,
Magnetic resonance imaging - MRI,
Pozytonowa tomografia emisyjna / Tomography - PET / CT,
X-RAY X-ray,
Digital subtraction angiography - DSA,
Diffusion spectrum imaging - DSI,
Image fusion,
Nuclear medicine,
Technical data:
Size: 30 ',
Color depth: 1.07 billion,
Resolution: 2560 x 1600,
Pixel size: 0.2505 x 0.2505 mm,
Brightness: 700 cd / m2,
Contrast: 1000: 1,
Viewing angles: 178 ° pion / 178 ° level,
Response time: 6 ms,
Signal inputs:
DisplayPort: 20 pins * 2
DVI-D: 24 pins * 2
VGA: D-Sub 15 pins
With table base,
Dimensions: 699 (W) x 513 ÷ 633 (H) x 251 (Depth) mm,
Own weight: 13 Kg,
Power supply: DC 24 V,
Power consumption: maximum 80 W,
Energy saving: EPA, VESA,
MTBF: > 50,000 hours,
Certificates: EN 60601, CCC, FCC, CE, UL,
Has a current Technical Passport issued,
36-month warranty

