



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



### **Short-wave diathermy Autotherm 390 (Reconditioned)**

Used, perfect technical condition,

U.S.A. production (2012).,

Power supply: 230 V,

Frequency: 50-60 Hz,

Autotherm 390 is a modern shortwave diathermy for high frequency electromagnetic field therapy. The device can work in a continuous or pulse mode, which allows you to obtain a thermal or thermal effect,

Advantages of Autotherm 390:

Simple service,

Readable control panel,

Continuous and impulse work,

Identification of connected electrodes,

Microprocessor controlled device with sound signaling,

Automatic reset of the output power after the end of the therapy,

Module for automatic camera testing when switching on,

Acoustic signal for ending therapy,

Technical data:

Generated frequency: 27.12 MHz

Type of emission: continuous, pulsating,

Power:

Continuous emission: 0 - 100 W,

Pulsation emission: 0 - 200 W,

Pulse frequency: 10 Hz, 20 Hz, 50 Hz, 100 Hz, 400 Hz,

Pulse duration: 65  $\mu$ s, 100 $\mu$ s, 200 $\mu$ s, 300 $\mu$ s, 400 $\mu$ s,

Clock, exposure time: 1-30 minutes,

Included:

2 rubber applicators with dimensions of 12 x 18 cm,

6 felt pads for rubber applicators,

2 cotton covers,

2 elastic tapes measuring 122 cm for fixing rubber electrodes,

Dimensions: 36 cm (width) x 42 cm (D) x 18 cm (height),

Weight: 6.8 kg,

User manual in Polish,

Has a current Technical Passport issued,

1-month warranty

