



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybierski
+48 882 762 006
kw@praiston.pl



Critical Parameter Analyzer ABL800 (New)

- New,
- Manufactured in Denmark,
- The FLEXQ module automatically identifies, mixes and automatically measures three samples simultaneously, so automatic identification and mixing reduces the risk of user errors,
- Measures pH, pO₂, pCO₂, Na⁺, K⁺, Ca⁺⁺, Cl⁻, glucose, lactate, nbils, tHb, FO₂Hb, FCO_{Hb}, FHHb, FMetHb,
- Time to obtain the first result 1.33 min,
- Reagent modules,
- Touchscreen display,
- Internal printer,
- Built-in sequential video system,
- Automatic sample medium recognition,
- Automatic calibration,
- Possibility of cooperation with laboratory information system,
- Barcode reader,
- Windows XP operating system,
- Automatic result processing, test ordering and data transmission to LIS,
- Automatic and customized data recording,
- Backup,
- CD-ROM drive for reading/writing,
- USB port,
- Dimensions: 55 x 53 x 70 (cm),
- Weight: 35kg,
- Valid technical passport issued,
- Warranty 6 months,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

