



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Diabetic analyzer Care Diagnostica EcoMatic (Reconditioned)

Used, very good technical condition,

German production (2013),

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 30 VA,

The analyzer is used to determine the concentration of glucose or lactate in whole blood,

Measurement method: enzymatic-amperometric,

Short measurement time and camera performance up to 130 tests per hour ensure quick results,

Touch, color LCD screen,

Built-in thermal printer,

Sample detection: automatic,

Plate capacity: 50 samples, 1 item 'cito', 1 standard, 3 checks,

Calibrations: start, cyclical, periodic,

Built-in connectors: RS232, USB,

Selection of units: mg / dL or mmol / L,

Memory: up to 1000 results (with date and time of examination),

Software: in Polish, English or German,

The volume of the sample: 20 microliters (whole blood, plasma, serum),

Sample collection: end-to-end capillary tube + 1.5 / 2.0 ml tube, with reagent,

Measuring range: Glucose 9 - 900 mg / dl,

Precision: <3% (at 12 mmol / liter); systematic error <= 2.5%,

Weight: 6 kg,

Overall dimensions: 40 x 43 x 24.5 cm,

Has a current Technical Passport issued,

1-month warranty

