



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Urine Analyzer AMP USR-500 (Reconditioned)

Used, very good technical condition,

Austrian production (2011),

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 40 W,

AMP USR-500 urine strip reader is a functional and efficient analyzer for automatic reading of AMP urine stripes in laboratories with high sample throughput up to 514 strips per hour,

Operation via an intuitive touch screen,

It has a built-in thermal printer,

Software in English,

Automatic calibration,

Automatic transport of reading strips,

Automatic removal of used strips to the built-in container,

Compensation of intense urine,

Opportunity to cooperate with a laboratory information system,

Determination of up to 11 parameters: URO, BIL, KET, BLD, PRO, NIT, LEU, GLU, SG, pH, ASA,

Type of test strips: AMP Urine Strips US-8, AMP Urine Strips US-10, AMP Urine Strips US-11,

Measurement method: colorimetric with the use of reflected light detection,

Wavelengths: 525, 572, 610 and 660 nm,

Memory of results: 2000 determinations + 50 quality control results,

Communication: RS-232 serial interface,

Included:

Barcode reader,

Weight: 7.5 kg,

Dimensions: 305 x 380 x 378 mm,

Has a current Technical Passport issued,

1-month warranty,

Cat 02

